

Software Development Kit

To allow users of the [DSX Suite](#) to expand the suite with their own code, a DSX Software Development Kit, the (**DSX SDK**), is included as part of the suite. The DSX SDK comes in its own installer file and needs to be installed separately from the DSX Suite. Users can add custom algorithms to the DSX Suite using the DSX SDK to compile and build additional dynamic library files (dll files) that most DSX applications will load during start-up.

Currently the SDK enables our DSX users to develop their own LCS algorithms that can be applied in [Orient3D](#). We plan to expand the SDK in the future.

DSX SDK Software Requirements

In order for the generated dll files to work properly with the DSX applications they need to be compiled and build with the same compiler that is used for DSX Suite. Currently this is: **Microsoft Visual Studio 2019**.

Debugging the newly developed code requires the **DSX Suite** needs to be installed as well. The use of the Qt and/or VTK libraries is not required, but if they are used, their versions need to match the ones with which the DSX Suite is developed.

Current **Qt** version: 5.13.1

Current **VTK** version: 8.2

Downloading the DSX SDK

HAS-Motion software is all downloaded over the Internet at service.has-motion.com. Customers are provided login information to access the download page as part of their initial email.

After logging into your customer account, you will see a screen similar to this:

C-Motion Research Branches

Products + Support + About + Free Downloads Download My Profile

Software Downloads | My Profile | Activated Licenses | Lab Contacts | Documentation | Logout

Software Downloads

Software Purchased	License Allowed	License Activated	Support Expires	Files to Download
Orient3D Professional	10	10	Never	Orient3D_Devel_v01_v01.01.01.exe Orient3D_Devel_v01_v01.01.01.exe
PointCloud	1	0	Never	PointCloud_SDK_Setup.exe License Key: 9514228A886A8877
Fused3D	2	2	Never	Fused3D_SDK_Setup.exe License Key: 8948249321373023
CellTracer Application	4	0	Never	CellTracer_Devel_v01_v01.01.01.exe CellTracer_Devel_v02_v01.01.01.exe
SDK Evaluation	1	0	2011-10-31	C-Motion SDK Evaluation_Setup.exe License Key: 712308801588-0196

Plug-in Modules

Module	Download
RT Server	RTServer_V00001_Setup.exe
Qualtrics	Qualtrics.dll

Activated Licenses

Product	License Key	File Name	Ever Activated	Comments
---------	-------------	-----------	----------------	----------

Click the **DOWNLOAD** Buttons to save the DSX SDK installation file to your computer. Most web browsers will offer you a choice of where to save the downloaded file; make sure you remember

where you save the downloaded file. Run the installation executable file you downloaded.

We recommend saving at least one prior installation file of the DSX SDK that works for you, as older installation files may be removed from the HAS-Motion download page without notification.

Installing the DSX SDK

When you run the executable file that you downloaded you will be guided through the installation steps. Your computer may ask you to confirm that you trust the source of the program you are installing.

The installer will install

1. a number of header files in the *include* folder
2. one source file (CMGlobalTime.cpp) in the *src* folder
3. 4 static library files in the *lib* folder
4. a working example in the *Example* folder.

It will also add the **CM_DSXSDK_DIR** variable to the computer's registry. This environment variable enables the DSX applications to find your dll files. This will be a system-wide environment variable if the installer is run with administrator privileges, otherwise it will only be available to the current user.

When the Completed dialog appears the DSX SDK has been successfully installed on your computer.

Using the DSX SDK

The easiest way to get started is to copy the contents of the *Example* folder in the SDK installation folder to a location where you have write permission. The names of the three *.vcxproj files can be changed in Windows Explorer (e.g. into MyDll.vcxproj, MyDll.vcxproj.filters, and MyDll.vcxproj.user) as long as all three files have the same name and their extensions remain unchanged. The project can be loaded into Visual Studio by double clicking on the *.vcxproj file or by starting Visual Studio, selecting 'Open a project or solution' and selecting the .vcxproj file in the file dialog.

The names of the header and source files can be changed in Visual Studio. New files (with new C++ classes) can be added as well. On exit Visual Studio will give the opportunity to save the *.sln (solution) file. This is not required, but is a means of keeping multiple VS projects (i.e., multiple dlls) together.

Export directive

The LCSExampleModule.h (or MyDllModule.h) file defines a macro that defines the export directive. When the name of the macro is changed, the LCSExampleTransform.h (or MyTransform.h) and LCSExampleFactory.cpp (or MyFactory.cpp) files need to be edited accordingly as well. Make sure the MYLCS_EXPORTS entry matches the preprocessor definition in the properties of the project for both Debug and Release mode.

Algorithm class

To implement a new LCS algorithm, one would create a C++ class like the LCSExampleTransform class; the new class would inherit from CMLCSTransformBase, and its constructor and internalUpdate() are the only functions that need to be overloaded.

- The *Description* variable is the identifier that is displayed in the O3D gui and is set using the setDescription() function. Internally, this string is used for the instantiation of the correct algorithm class object.
- The *m_Side* variable enables the implementation of both the proximal and distal LCS algorithms for a particular bone in one C++ class.
- The *m_Landmarks* variable stores the mandatory and optional landmarks used by the algorithm. Each landmarks constructor takes a string for the name that is displayed in the application gui. By default a landmark is not optional; use the setOptional() function to change this.

The implementation of the algorithm goes into the internalUpdate() function. The function will use the object's surface model data (*m_Mesh*; for the CM3DFmtBase class see the CM3DFmt.h file in the \include\CMLib\CMCore\CM3DFmt folder in the DSX SDK installation folder) and the landmark data (*m_Landmarks*) to determine the elements of the 4×4 row-major **Matrix** member variable. The transformation matrix represents the transform from the object's CT coordinate system to its local coordinate system ([LCS](#)).

Factory

Additionally a factory class like the LCSExampleFactory is needed. The class is derived from cmFactoryBase and is used in the DSX applications to create the correct instantiation of a LCS algorithm class object from the description string that is chosen in the application. Unless the name of the factory class is changed, the header file does not need to be edited. The extern functions in the MyFactory.cpp are used by the DSX applications for some checks and the creation of the factory object during runtime. Make sure the correct class object gets created in the cmLoadFactory() function:

```
'%' extern "C" %'***'%MyDLL_EXPORT MyFactory%'***'%* cmLoadFactory()
{ return %'***'%MyFactory%'***'%::createInstance(); } %'
```

The CM_CREATE_CREATE_FUNCTION macro (defined in cmFactoryBase.h) is used to create the static functions that the factory uses to create the algorithm class objects and is used as follows:

```
'%' CM_CREATE_CREATE_FUNCTION(%'***'%MyTransform%'***'%) %'
```

In the constructor of the factory class the various algorithm classes are registered:

```
'%' %'***'%MyFactory%'***'%::%:'***'%MyFactory%'***'% (const
char* identifier) : cmFactoryBase(identifier)%'
'%'{ %'
```

```
'%' %'
```

```
%//% register built-in classes
```

```
**MyTransform** lcs_algo = **MyTransform**::createInstance();  
this->registerClass("CMLCSTransformBase", %%//%% changes for different  
algorithm types  
lcs_algo->getClassName(),  
lcs_algo->getDescription(),  
0,  
cmObjectFactoryCreate**MyTransform**);  
lcs_algo->deleteInstance();  
  
' '%%} %%'
```

NOTE: Depending on the situation, each factory uses either the string returned with *CMObjectbase::getClassName()* or *CMObjectbase::getDescription()* to create one of its registered class objects. In case of duplicates the first instance that is found is used.

NOTE: Ignore the warning C4275: non dll-interface class 'A' used as base for dll-interface class 'B' that Visual Studio generates.


Debugging your code

Debugging code created with the DSX SDK requires the DSX Suite to be installed as well. Because the applications in the DSX suite are release builds, you can only debug the release build of your dll. The *.vcxproj.user file in the DSX SDK Example installation folder ensures the Orient3D application is started when debugging is started (F5 key) in Release mode.

NOTE: Realize that debugging release code means certain parts of the code might have been removed from the optimized code. Some additional information can be found [here](#) and [here](#).

Including your dll in DSX applications

The DSX applications need to know in which folder to look for your dll file. To accomplish this create an environment variable **CM_DSX_USER_DIR** that points to this folder.



From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=other:dsx:software_development_kit

Last update: **2025/01/17 18:43**

