

Comments can be added to .cen files with the 'Edit Comment' button before or after the recording.

[Capture_Browse.jpg](#)

The comment is stored within the cen file, so each recording can be annotated individually.

[editcomment.png](#)

The comment can be added prior to the recording, or immediately after the recording.

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

<https://has-motion.com/wiki/doku.php?id=other:capture2d:comments&rev=1718647979>

Last update: **2024/06/17 18:12**

