

Use the browse button to select a good location in which to store your cen files. The default path is the Capture application folder, which isn't writable on Windows 7 and isn't a smart place to store data (fix coming soon). Also be sure to give your files good names (i.e., cal.cen for calibration trial, jump.cen for a trial that involves jumping, etc.).

Retrieved from ""

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:capture2d:where_to_save_on_disk&rev=1718643908

Last update: **2024/06/17 17:05**

