

use the browse button to select a good location in which to store your cen files. the default path is the capture application folder, which isn't writable on windows 7 and isn't a smart place to store data (fix coming soon). also be sure to give your files good names (i.e., cal.cen for calibration trial, jump.cen for a trial that involves jumping, etc.).

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=other:capture2d:where_to_save_on_disk&rev=1718801144

Last update: **2024/06/19 12:45**

