

# Convert3D\_Overview

## Overview

Convert3D allows users to convert files from one format to another. The file types that can be converted are:

- C3D to CMO
- TXT to V3D TXT

Value to users:

- Allows users to create .cmo files for Inspect3D without using Visual3D
- Allows users to create .txt files that can be read by Visual3D

## Wiki Navigation

- [Installation](#)
- [Documentation](#)
- Tutorials
- Troubleshooting

From:  
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=other:convert3d:convert3d\\_overview&rev=1720792519](https://has-motion.com/wiki/doku.php?id=other:convert3d:convert3d_overview&rev=1720792519)

Last update: **2024/07/12 13:55**

