

Convert3D_Overview

Overview

Convert3D allows users to convert files from one format to another. The file types that can be converted are:

- C3D to CMO
- TXT to V3D TXT

Value to users:

- Allows users to create .cmo files for Inspect3D without using Visual3D
- Allows users to create .txt files that can be read by Visual3D

Wiki Navigation

- [Installation](#)
- [Documentation](#)
- Tutorials
- Troubleshooting

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:convert3d:convert3d_overview&rev=1721149283

Last update: **2024/07/16 17:01**

