


Locate3D Settings

The Located3D settings dialog allows users to view and modify application settings.



- **Euler** the order of Euler rotations in the body pose maps. XYZ is the default, but in some situations it may be helpful to change the order to avoid gimbal lock.
- **Color Theme** you can choose between a light or a dark color theme for the GUI.

See Also

Go back to the [Locate3D Overview](#).

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=other:dsx:locate3d:locate3d_settings

Last update: **2025/05/05 21:03**

