


# Object Configuration

Orient3D's Object Configuration dialog allows users to select which subject object (either a [surface file](#) or [ROI surface file](#)) is displayed in the [3D view](#).



|          |                                                                                                                                                            |
|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Subject: | 09                                                                                                                                                         |
| Session: | 2017_04_11                                                                                                                                                 |
| Object:  | <input type="checkbox"/> left_femur<br><input type="checkbox"/> left_tibia<br><input type="checkbox"/> right_femur<br><input type="checkbox"/> right_tibia |

- **Subject** This field is automatically populated with the ID of the subject when you load a subject file.
- **Session** This drop down menu is automatically populated with all of the sessions available for the current subject.
- **Object** This list is automatically populated with all of the subject objects when you select a session. If the object has a *Surface File* or *ROI Surface File* specified for it, the object is enabled in the list so you can select it. The object that is selected is displayed in the 3D view. You can select only one object at a time.

## See Also

Go back to the [Orient3D Overview](#).

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=other:dsx:orient3d:object\\_configuration](https://has-motion.com/wiki/doku.php?id=other:dsx:orient3d:object_configuration)

Last update: 2025/05/07 13:15

