


Object Configuration

A screenshot of a software window titled "Object Configuration". The window has a light gray background and a purple title bar. It contains three fields: "Subject:" with the value "09", "Session:" with a dropdown menu showing "2017_04_11", and "Object:" with a list of four items: "left_femur", "left_tibia", "right_femur", and "right_tibia", each preceded by an unchecked checkbox.

Subject:	09
Session:	2017_04_11
Object:	☐ left_femur ☐ left_tibia ☐ right_femur ☐ right_tibia

- **Subject** This field is automatically populated with the ID of the subject when you load a subject file.
- **Session** This drop down menu is automatically populated with all of the sessions available for the current subject.
- **Object** This list is automatically populated with all of the subject objects when you select a session. If the object has a *Surface File* or *ROI Surface File* specified for it, the object is enabled in the list so you can select it. The object that is selected is displayed in the 3D view. You can select only one object at a time.

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

https://has-motion.com/wiki/doku.php?id=other:dsx:orient3d:object_configuration&rev=1721231057

Last update: 2024/07/17 15:44

