


Orient3D Settings

The Orient3D Settings dialog allows the user to modify application settings.



- **Landmark Scale** Changes the size of the landmark spheres and text labels in the 3D window.
- **Color Theme** You can choose between a light and dark color theme.

See Also

Go back to the [Orient3D Overview](#).

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=other:dsx:orient3d:orient3d_settings

Last update: **2025/05/07 18:24**

