


PlanDSX Settings

The PlanDSX Setting dialog allows the user to change two application-wide default values.



- Units: The units for translations and distances in the scene. **Note:** This is not implemented.
- Color Theme: You can choose between a light or a dark colour theme for the GUI.

[Back to PlanDSX Overview](#)

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:dsx:plandsx:plandsx_settings

Last update: **2025/01/10 21:19**

