

Surface3D_Settings



- **Undo Limit** The maximum number of label editing actions that can be undone. The default value is 10.
- **Paintbrush Radius** The radius (in voxels) of the circular paintbrush that is used for painting and erasing segmentation labels. The default value is 5. The radius can also be increased or decreased using the mouse wheel, after pressing **p** to paint or **e** to erase.
- **Text Scale** The scale factor for POI and landmark text labels. The default value is 0.75.
- **Color Theme** You can choose between a light or a dark color theme.

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:dsx:surface3d:surface3d_settings&rev=1720792559

Last update: **2024/07/12 13:55**

