

# DSX\_Configuration



- **Subject** This field is automatically populated with the ID of the subject when you load a subject file.
- **Session** This drop down menu is automatically populated with all of the sessions available for the current subject.
- **Configuration** This drop down menu is automatically populated with all of the configurations available for the currently selected session.
- **Trial** This drop down menu is automatically populated with all of the trials available for the currently selected configuration.
- **Object** This list is automatically populated with all of the objects in the currently selected trial. Each object is enabled if it is trackable, meaning that it has a 3D image file associated with it, for generating DRRs.

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

[https://has-motion.com/wiki/doku.php?id=other:dsx:x4d:dsx\\_configuration&rev=1721149309](https://has-motion.com/wiki/doku.php?id=other:dsx:x4d:dsx_configuration&rev=1721149309)

Last update: **2024/07/16 17:01**

