


DSX_Configuration



DSX Configuration

Subject: KneeSubject

Session: Day1

Configuration: Left Level Walk

Trial: LeftStand

Object:

- ☒ All
- ☒ left_femur
- ☒ left_tibia

- **Subject** This field is automatically populated with the ID of the subject when you load a subject file.
- **Session** This drop down menu is automatically populated with all of the sessions available for the current subject.
- **Configuration** This drop down menu is automatically populated with all of the configurations available for the currently selected session.
- **Trial** This drop down menu is automatically populated with all of the trials available for the currently selected configuration.
- **Object** This list is automatically populated with all of the objects in the currently selected trial. Each object is enabled if it is trackable, meaning that it has a 3D image file associated with it, for generating DRRs.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:dsx:x4d:dsx_configuration&rev=1721228978

Last update: **2024/07/17 15:09**

