

Settings



- **CPU Threading** With two CUDA-capable graphics cards, object tracking can be made faster when each card is controlled by its own CPU thread. To assist with customer support, C-Motion staff might occasionally ask to disable CPU threading when debugging problems with graphics cards.
- **X-Ray Opacity** The opacity of the X-ray images in the 2D and 3D windows. Their opacity should be less than 1.0 so that the DRRs can be seen superimposed on them. The default value is 0.5.
- **Simulated Annealing Temp** The initial temperature for the ASA simulated annealing algorithm. The higher the initial temperature, the more slowly the search range will narrow during the optimization. The default value is 2000.
- **Color Theme** You can choose between a light and dark color theme for the graphical interface.

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
<https://has-motion.com/wiki/doku.php?id=other:dsx:x4d:settings&rev=1720792575>

Last update: **2024/07/12 13:56**

