

Line_Styles

Inspect3D Line Styles

Specify line style you want to define. You can specify a line style for:

- Query (Group)
- Workspace
- Condition
- Exclusion
- Classes
- Vertical Line

Line Style

- Thickness: thickness of line
- Line Style: style of line
- Scatter Style: highlight each data point (ex. "X")
 - Scatter Size: size of data point

NOTE: You can draw just the line, or just the data points by setting Line Style or Scatter Style to None.

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:inspect3d:documentation:definitions:line_styles&rev=1720792593

Last update: **2024/07/12 13:56**

