

Line Styles

Inspect3D Line Styles

Specify line style you want to define. You can specify a line style for:

- Query (Group)
- Workspace
- Condition
- Exclusion
- Classes
- Vertical Line

Line Style

- **Thickness:** thickness of line
- **Line Style:** style of line
- **Scatter Style:** highlight each data point (ex. "X")
 - **Scatter Size:** size of data point

NOTE: You can draw just the line, or just the data points by setting Line Style or Scatter Style to None.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:inspect3d:documentation:definitions:line_styles&rev=1734456427

Last update: **2024/12/17 17:27**

