


Report3D_CMO_Library_Browser

Introduction

The purpose of the CMO Library Browser Window is to show the user a summary of the CMO files, c3d files, Tags and Events that are available to them.

CMO Library Browser Parts



Loaded CMO Files - This drop down menu has a list of all the CMO files that you have loaded into the workspace. The CMO file that is selected has the C3D files that it contains listed in the C3D Files box.

C3D Files - This box lists the C3D Files which are contained in the CMO file selected in the Loaded CMO Files drop down menu.

Library Tags - This drop down menu has a list of all the Tags in the CMO file selected in the Loaded CMO Files drop down menu.

Library Events - This drop down menu has a list of all the Events in the CMO file selected in the Loaded CMO Files drop down menu.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:report3d:documentation:report3d_cmo_library_browser&rev=1721228993

Last update: 2024/07/17 15:09

