

Report3D CMO Library Browser

Introduction

The purpose of the CMO Library Browser Window is to show the user a summary of the CMO files, c3d files, Tags and Events that are available to them.

CMO Library Browser Parts

CMO Library Browser

Loaded CMO Files

s043

C3D Files

functHipL

functHipR

fwalk0001

fwalk0002

fwalk0003

Library Tags

FAST

Library Events

LHS

Loaded CMO Files - This drop down menu has a list of all the CMO files that you have loaded into the workspace. The CMO file that is selected has the C3D files that it contains listed in the C3D Files box.

C3D Files - This box lists the C3D Files which are contained in the CMO file selected in the Loaded CMO Files drop down menu.

Library Tags - This drop down menu has a list of all the Tags in the CMO file selected in the Loaded CMO Files drop down menu.

Library Events - This drop down menu has a list of all the Events in the CMO file selected in the Loaded CMO Files drop down menu.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:report3d:documentation:report3d_cmo_library_browser&rev=1721229724

Last update: 2024/07/17 15:22

