

Report3D_Page_Previews

Page Previews



The Page Previews window allows the user to see thumbnails of the pages in the report. The page that is active in the view area is highlighted in blue in the Page Preview window. The window can be closed using the X in the top right hand corner. It can also be removed from the Report3D screen and viewed separately by either clicking on the double squares in the top right hand corner or by dragging the window using the bar at the top.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:report3d:documentation:report3d_page_previews&rev=1720790837

Last update: 2024/07/12 13:27

