

Report3D Page Previews

Page Previews

Page Previews

Page ID: {d12b736b-0d22-4087-8d2f-54659acdb5cd}

Add

Remove

The Page Previews window allows the user to see thumbnails of the pages in the report. The page that is active in the view area is highlighted in blue in the Page Preview window. The window can be closed using the X in the top right hand corner. It can also be removed from the Report3D screen and viewed separately by either clicking on the double squares in the top right hand corner or by dragging the window using the bar at the top.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:report3d:documentation:report3d_page_previews&rev=1721230915

Last update: 2024/07/17 15:41

