

report3d wiki overview page]]. ===== insert object ===== r3daddtextitem.png ****add text item**** - this button adds a text box to the current page of the report. r3daddtableitem.png ****add table item**** - this button adds a table to the current page of the report. r3dadd2dgraphitem.png ****add 2d graph item**** - this button adds a 2d graph to the current page of the report. r3daddpictureitem.png ****add picture item**** - this button adds a picture to the current page of the report. r3daddmovieitem.png ****add movie item**** - this button adds a movie to the current page of the report. r3daddrectangleitem.png ****add rectangle item**** - this button adds a rectangle to the current page of the report. r3daddellipseitem.png ****add ellipse item**** - this button adds an ellipse to the current page of the report. ===== move objects ===== r3dbringitemtofront.png ****bring item to front**** - this button moves the selected item to the front. r3dbringitemforward.png ****bring item forward**** - this button moves the selected item forward. r3dsenditembackward.png ****send item backwards**** - this button moves the selected item backwards. r3dsenditemtoback.png ****send item to back**** - this button moves the selected item to the back.

} } } }

From: <https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:report3d:documentation:report3d_tool_bar&rev=1718804100

Last update: **2024/06/19 13:35**

