

open the file drop down menu at the top left hand corner of the screen.\\ \\ select export as pdf.\\ \\ you will be asked where you would like to save the pdf and what you would like to call it.\\ \\ \\ |epstep2.png| ===== step 3. verify pdf ===== open the pdf document that you have saved and verify that the content is there. each of the report pages seen in the center of the screen will have been transferred to the pdf document.

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=other:report3d:report3d_export_pdf&rev=1718804096

Last update: **2024/06/19 13:34**

