

Report3D_Export_PDF

Purpose

WARNING: The Report3D product was never finished, so the title is wrong and the product is not available.

However, if you want to simply print a Visual3D report as a PDF, you can do this:

- 1) If you have Adobe Acrobat Pro, jump to step 2. Otherwise, download and install any shareware PDF creation package, like the CutePDF Writer (www.cutepdf.com) and install it.
- 2) Open your Visual3D report and click the printer icon in the toolbar. Select the PDF tool as the output device. Give the PDF file a meaningful name.

Steps

Step 1. Open Report


Make sure that you have a report open. You can make it from scratch, alter an existing report, or load an existing report.

Step 2. Export and Save

Open the file drop down menu at the top left hand corner of the screen.

Select Export as PDF.

You will be asked where you would like to save the pdf and what you would like to call it.

The screenshot shows the 'Report3D v1.0.0.204 - IORgaitFoot_Demo.r3d' application window. The 'File' menu is open, displaying options: New (Ctrl+N), Open... (Ctrl+O), Load CMO Library..., Save (Ctrl+S), Save As..., Export as PDF... (highlighted), Print... (Ctrl+P), and Exit. The application's toolbar is visible on the right, featuring icons for opening, saving, and printing, along with a text field containing 'Sub_003'.

Step 3. Verify PDF

Open the PDF document that you have saved and verify that the content is there. Each of the report pages seen in the center of the screen will have been transferred to the PDF document.

Last
update:
2024/07/17 15:09 other:report3d:report3d_export_pdf https://has-motion.com/wiki/doku.php?id=other:report3d:report3d_export_pdf&rev=1721228993

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:report3d:report3d_export_pdf&rev=1721228993

Last update: **2024/07/17 15:09**

