

Report3D Export PDF

Purpose

WARNING: The Report3D product was never finished, so the title is wrong and the product is not available.

However, if you want to simply print a Visual3D report as a PDF, you can do this:

- 1) If you have Adobe Acrobat Pro, jump to step 2. Otherwise, download and install any shareware PDF creation package, like the CutePDF Writer (www.cutepdf.com) and install it.
- 2) Open your Visual3D report and click the printer icon in the toolbar. Select the PDF tool as the output device. Give the PDF file a meaningful name.

Steps

Step 1. Open Report


Make sure that you have a report open. You can make it from scratch, alter an existing report, or load an existing report.

Step 2. Export and Save

Open the file drop down menu at the top left hand corner of the screen.

Select Export as PDF.

You will be asked where you would like to save the pdf and what you would like to call it.



Step 3. Verify PDF

Open the PDF document that you have saved and verify that the content is there. Each of the report pages seen in the center of the screen will have been transferred to the PDF document.

Last
update:
2024/07/17 15:44 other:report3d:report3d_export_pdf https://has-motion.com/wiki/doku.php?id=other:report3d:report3d_export_pdf&rev=1721231078

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=other:report3d:report3d_export_pdf&rev=1721231078

Last update: **2024/07/17 15:44**

