

Visual3D_Server_Plugins_Overview

AMTI Gen5

The AMTI Gen5 plugin supports the [Gen5](#) USB Driver 1.0 (June 2010) and newer.

AMTI Gen5 Type6

The AMTI Gen5 Type6 plugin supports the same hardware/software as the AMTI Gen5 plugin but allows force plates to be formatted as a Type 6 force plate.

Bertec AM6X00

The Bertec AM6X00 plugin supports [Bertec Force Plates](#) with software versions 1.0.0 or newer.

C3D File

The C3D file plugin supports c3d files formatted according to [c3d.org](#).

CEN File

The CEN File plugin supports .cen files which were recorded by 2D cameras and which have the position, width and height of centroids recorded.

MAC_Cortex

The MAC Cortex plugin supports [Cortex v3.x and v4.1 or newer](#).

NaturalPoint Motive

Myon

The Myon plugin is under development for [Myon products](#).

NDI First Principles

The NDI First Principles plugin supports [First Principles real-time protocol v1.2.2 or newer](#).

NP OptiTrack

The NP OptiTrack plugin provides support for the following [NaturalPoint™ OptiTrack™](#) camera systems:

Cameras

- [Flex 100:R1™](#), [Flex 100:R2™](#), [Flex 100:R2 FS™](#)
- [Flex 13™](#)
- [S250e™](#)

Hubs

- [OptiHub™](#), [OptiHub 2™](#)
- [eSync™](#)

This plugin is up to date with the camera SDK version 1.7.5.

Qualisys QTM

The Qualisys QTM plugin supports [QTM real-time protocol](#) v1.8 or newer.

Vicon DataStream

The Vicon DataStream plugin supports [Nexus v1.7.x, v1.8.x, 2.0, and 2.1](#) and [Tracker 1.3, and 2.0](#).

x-io

The x-io plugin supports [x-IMUs](#).

From: <https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link: https://has-motion.com/wiki/doku.php?id=other:visual3dserver:plugins:visual3d_server_plugins_overview&rev=1720792676

Last update: **2024/07/12 13:57**

