

By right-clicking on a plot, users can select "Legend" from the context menu. The Legend menu will enable options once the user selects "Legend Visible". The legend automatically populates the signal names based on the group and if it is the raw data or the dispersion. By making the legend visible, users can select the alignment of the legend (in any corner of the plot).



From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

[https://has-motion.com/wiki/doku.php?id=sift:application:plot\\_legend\\_dialog&rev=1721157827](https://has-motion.com/wiki/doku.php?id=sift:application:plot_legend_dialog&rev=1721157827)

Last update: 2024/07/16 19:23

