


Plot Legend Dialog

Right-clicking on a plot uses a context menu to “Legend” from the context menu. The Legend menu will enable options once the user selects “Legend Visible”. The legend automatically populates the signal names based on the group and if it is the raw data or the dispersion. By making the legend visible, users can select the alignment of the legend (in any corner of the plot).



Last
update:
2024/07/17 15:10 sift:application:plot_legend_dialog https://has-motion.com/wiki/doku.php?id=sift:application:plot_legend_dialog&rev=1721229005

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=sift:application:plot_legend_dialog&rev=1721229005

Last update: **2024/07/17 15:10**

