


Plot Legend Dialog

Plot legends can be added to any plot on the explore page and any of the four plot types: Signal-Time, Metric, Signal-Signal, and Metric-Metric. They can also be added to most plots on the Analyse page, including PCA and SPM plots.



By right clicking on a plot, users can select “Legend” from the context menu. The Legend menu will enable options once the user selects “Legend Visible”. The legend automatically populates the signal names based on the current [Display Style](#), as well as the plotted type (workspace mean, dispersion, etc.). By making the legend visible, users can select the alignment of the legend (in any corner of the plot).

☒ Legend Visible

Legend Alignment

Top Right

ANKLE_ANGLE_X Workspace Means

ANKLE_ANGLE_X Workspace Std Devs

Cancel

Save

☒ Legend Visible

Legend Alignment

Top Right

Legend Size

Small

AnkleAngleX_::Sub05\Sub05_Workspace.c

◀

▶

Cancel

Save

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=sift:application:plot_legend_dialog&rev=1731681719

Last update: 2024/11/15 14:41

