


Select Font Dialog

The Select Font Dialog is a font style editor for plot titles and axes throughout Sift. It is generally accessed through the context menu on a graph, but is also available through the [Plot Options Dialog](#)



- **Font:** The currently selected font, and all selection box containing all available fonts
- **Font Style:** The currently selected font styles, and all available font styles, **Bold**, *italic* etc.
- **Size:** The currently selected font size, and all available font sizes ranging from 6 to 72
- **Effects:** Contains two check-boxes to enable strike through (Example) and underline (Example) or both (Example).
- **Writing System:** A selection of different character sets.
- **Sample:** An example box previewing what the selected font will look like.
- **OK:** Close the dialog and apply all font choices
- **Cancel:** Close the dialog and do not apply font choices

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=sift:application:select_font_dialog

Last update: **2024/11/15 16:25**

