

| the select font dialog is a font style editor for plot titles and axes throughout sift.

sift_font_select.png\\ * **font**: all selection box containing all available fonts * **font style**: all available font styles, **bold**, *italic* etc. * **size**: all available text sizes ranging from 6 to 72 * **effects**: contains two check-boxes to enable strike through (example) and underline (example) or both (example). * **writing system**: a selection of different character sets. * **sample**: an example box previewing what the selected font will look like. * **ok**: close the dialog and apply all font choices * **cancel**: close the dialog and do not apply all font choices

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=sift:application:select_font_dialog&rev=1718804107

Last update: **2024/06/19 13:35**

