




- **Font:** All selection box containing all available fonts
- **Font Style:** All available font styles, **Bold**, *italic* etc.
- **Size:** All available text sizes ranging from 6 to 72
- **Effects:** Contains two check-boxes to enable strike through (Example) and underline (Example) or both (Example).
- **Writing System:** A selection of different character sets.
- **Sample:** An example box previewing what the selected font will look like.
- **OK:** Close the dialog and apply all font choices
- **Cancel:** Close the dialog and do not apply all font choices

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

https://has-motion.com/wiki/doku.php?id=sift:application:select_font_dialog&rev=1721229734

Last update: 2024/07/17 15:22

