


# Select Font Dialog

The Select Font Dialog is a font style editor for plot titles and axes throughout Sift. It is generally accessed through the context menu on a graph, but is also available through the [Plot Options Dialog](#)



- **Font:** All selection box containing all available fonts
- **Font Style:** All available font styles, **Bold**, *italic* etc.
- **Size:** All available text sizes ranging from 6 to 72
- **Effects:** Contains two check-boxes to enable strike through (Example) and underline (Example) or both (Example).
- **Writing System:** A selection of different character sets.
- **Sample:** An example box previewing what the selected font will look like.
- **OK:** Close the dialog and apply all font choices
- **Cancel:** Close the dialog and do not apply all font choices

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

[https://has-motion.com/wiki/doku.php?id=sift:application:select\\_font\\_dialog&rev=1731687824](https://has-motion.com/wiki/doku.php?id=sift:application:select_font_dialog&rev=1731687824)

Last update: 2024/11/15 16:23

