

tutorials page]]. if you have a problem that has not been covered below or in the sift tutorials feel free to e-mail the support team at support@has-motion.ca. \\ ===== working data ===== all grouped data ===== all data/workspaces that were queried, regardless of what is plotted. ===== all plotted groups ===== only perform work on the groups that are plotted on all plots. work can be performed on: * all workspaces for those groups (regardless on if they are plotted or not) * plotted workspaces (works on entire workspace, even if only one trace/c3d is plotted) ===== all workspaces ===== only perform work on the groups that are plotted in the active plot. * all workspaces for those groups (regardless on if they are plotted or not) * plotted workspaces (only works on the plotted data, so if only one trace/c3d is plotted, that is what will be used) ===== all highlighted data ===== only data that is highlighted ===== graph is missing data ===== if you plot data and there are empty data points, you may end up with a straight line at -9999999: datamissinggraph.png most calculations (mean, standard deviation, etc) will handle empty data points without an issue, so you can choose not to plot data not found: * in the [[sift:application:plot_options_dialog|options dialogue]] * uncheck **"show data not found"** * click **"refresh plot"** you can also choose to exclude this data by selecting this data, right clicking -> exclude -> exclude trace (raw data). adjusteddatamissinggraph.png ****note:**** most calculations handle data not found, however, the input to pca should not contain any missing data points. }}

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

<https://has-motion.com/wiki/doku.php?id=sift:documentation:troubleshooting&rev=1718804114>

Last update: **2024/06/19 13:35**

