

Free CMO Viewer

C-Motion provides a CMO viewer for users that do not have a Visual3D license.

The viewer allows you to open cmo files, review the model definition, see the original data, model based signals, any derived signals and their processing history.

The viewer allows you to see the Visual3D report.

In this way you can provide your biomechanical analyses and reports for others to see.

What the viewer does not do, is to allow you to edit the data.

The Free CMO viewer can be downloaded [Here\]](#)

A video of the installation process and an overview of the functionality can be found [Here\]](#)

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:cmo_viewer:free_cmo_viewer&rev=1721230922

Last update: **2024/07/17 15:42**

