

Subject Prefix: In the C3D file, there is a parameter to specify the subject prefix.

* This should be specified by your motion capture system.

* The subject prefixes can be different from C3D file to C3D file

When using multi-subject data, if the subject prefix is specified in the C3D file, you can select which subject to create segments for by specifying the subject prefix.

Retrieved from ""

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:definitions:subject_prefixes&rev=1718643946

Last update: **2024/06/17 17:05**

