

Subject Prefix: In the C3D file, there is a parameter to specify the subject prefix.	
* This should be specified by your motion capture system. * The subject prefixes can be different from C3D file to C3D file	
When using multi-subject data, if the subject prefix is specified in the C3D file, you can select which subject to create segments for by specifying the subject prefix.	

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:definitions:subject_prefixes&rev=1718805368

Last update: 2024/06/19 13:56

