

Subject_Prefixes

Subject Prefix: In the C3D file, there is a parameter to specify the subject prefix.

✖

* This should be specified by your motion capture system.
* The subject prefixes can be different from C3D file to C3D file

✖

When using multi-subject data, if the subject prefix is specified in the C3D file, you can select which subject to create segments for by specifying the subject prefix.

✖

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:definitions:subject_prefixes&rev=1719433304

Last update: 2024/06/26 20:21



Wiki Documentation Site - <https://has-motion.com/wiki/>