



# Subject\_Prefixes

**Subject Prefix:** In the C3D file, there is a parameter to specify the subject prefix.

\* This should be specified by your motion capture system.  
\* The subject prefixes can be different from C3D file to C3D file



When using multi-subject data, if the subject prefix is specified in the C3D file, you can select which subject to create segments for by specifying the subject prefix.



From:  
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:definitions:subject\\_prefixes&rev=1721229017](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:definitions:subject_prefixes&rev=1721229017)

Last update: **2024/07/17 15:10**

