

this command will open the the “all events” interface and “create events” interface if command is true.

if show_create_events_dialog = true, the command create_event_list can be populated to add additional events to the event list. similarly, clear_existing_event_list will remove chosen events.

open edit dialogs

```
! /show_all_events_dialog=false  
! /show_create_events_dialog=false  
! /create_event_list=  
! /clear_existing_event_list=false
```

example 1

open edit dialogs

```
/show_all_events_dialog=true  
! /show_create_events_dialog=false  
! /create_event_list=  
! /clear_existing_event_list=false
```



example 2

open edit dialogs

```
! /show_all_events_dialog=false  
/show_create_events_dialog=true  
/create_event_list=rhs+rtoe  
! /clear_existing_event_list=false
```



From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:dialogs:open_edit_dialogs&rev=1718804142

Last update: **2024/06/19 13:35**

