

**This command is deprecated. Instead export EMG signals from Noraxon directly to a C3D file and then import signals using the [Import\\_Signals\\_From\\_C3D\\_File](#) command.**

Imports data from [Noraxon EMG](#) files into an existing [C3D](#) file.

```
Import_Analog_Noraxon_EMG_Signals  
/FILE_NAMES=  
/ANALOG_RATE=  
;
```

[Back to Pipeline Commands Reference](#)

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:emg:import\\_analog\\_noraxon\\_emg\\_signals&rev=1718717259](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:emg:import_analog_noraxon_emg_signals&rev=1718717259)

Last update: **2024/06/18 13:27**

