

this command is deprecated. instead export emg signals from noraxon directly to a c3d file and then import signals using the [import_signals_from_c3d_file](#) command.

imports data from [noraxon emg](#) files into an existing [c3d](#) file.

```
import_analog_noraxon_emg_signals  
/file_names=  
/analog_rate=  
;
```

[back to pipeline commands reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:emg:import_analog_noraxon_emg_signals&rev=1718801256

Last update: **2024/06/19 12:47**

