

The default segment masses (relative to total body mass) are:

Segment Mass

foot	0.0145
shank	0.0465
thigh	0.10
hand	0.006
forearm	0.016
upper arm	0.028
pelvis	0.142
thorax/abdomen	0.355
head	0.081
trunk	0.497

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:modeling:segments:segment_mass&rev=1720792917

Last update: 2024/07/12 14:01

