

Segment Mass

Default Segment Mass

The default segment masses defined in Visual3D are based on [Dempster's regression equations](#).

The default segment masses (relative to total body mass) are:

foot	0.0145
shank	0.0465
thigh	0.10
hand	0.006
forearm	0.016
upper arm	0.028
pelvis	0.142
thorax/abdomen	0.355
head	0.081
trunk	0.497

Modifying Segment Mass

Segment Mass can be modified under the [segment properties](#) tab in Visual3D's model builder.

Segments | Segment Properties | IK Constraints | Right Thigh

Segment Name: Right Thigh

☒ Kinetic ☐ Kinetic But No Force Assignment ☐ Kinematic Only

Segment Mass: 0.1*Mass = 8.21 kgs

Segment Length: 0.420917 Meters

Segment Geometry: CONE

Inertial Values*		Center of Mass (Meters)*	
Expression	Value	Expression	Value
IXX:	0.12846	ML :	0*RTH_SEG_ 0
IYY:	0.12846	AP :	0*RTH_SEG_ 0
IZZ:	0.01735	Axial :	0.464927*RTI -0.195696

*May be modified for custom geometry

Model File: rthigh.v3g Browse

Material File: Browse

☐ Override Color Black

Rotate/Scale Graphic Model

Copy Override Color to Active Subject

Modify Segment Coordinate System

Apply

Build Model

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:modeling:segments:segment_mass&rev=1746214878

Last update: 2025/05/02 19:41

