

(modified from subject01_setup_scale.xml)

- the <scaletool> name has been changed to visual3d.
- the subject <mass> and <height> needs to be customized for each subject
- the <marker set> is no longer required and reference to this file has been removed
- <scaling_order> is set only to manualscale. (the manual scale factors for the segments will be written by visual3d to the visual3d_scale_scaleset.xml file.)
- the <measurementset> .xml file and its corresponding .trc file are no longer needed and reference to these files have been removed
- the <modelscaler> section has been removed since this is now done by visual3d
- the <markerplacer> section has been removed since this is now done by visual3d

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:opensim:setup_scale.xml&rev=1718804825

Last update: **2024/06/19 13:47**

