

Setup_Scale.xml

(modified from subject01_Setup_Scale.xml)

- The <ScaleTool> name has been changed to Visual3d.
- The subject <mass> and <height> needs to be customized for each subject
- The <Marker Set> is no longer required and reference to this file has been removed
- <scaling_order> is set only to manualScale. (The manual scale factors for the segments will be written by Visual3d to the Visual3d_Scale_ScaleSet.xml file.)
- The <MeasurementSet> .xml file and its corresponding .trc file are no longer needed and reference to these files have been removed
- The <ModelScaler> section has been removed since this is now done by Visual3d
- The <MarkerPlacer> section has been removed since this is now done by Visual3d

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:opensim:setup_scale.xml&rev=1721158062

Last update: **2024/07/16 19:27**

