

Event\_Import\_From\_Orthotrak

# Event\_Import\_From\_Orthotrak

To run this pipeline command:

Both the .c3d and Orthotrak generated .xls files must have the same name (except for the extension of course) and must be located in the same directory. Once the .c3d file has been loaded in Visual3D's workspace, the Event\_Import\_From\_Orthotrak can be executed from the pipeline. The events will automatically be imported in Visual3D and merged with its corresponding .c3d file.

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:event\\_commands:event\\_import\\_from\\_orthotrak&rev=1719433479](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:event_commands:event_import_from_orthotrak&rev=1719433479)

Last update: 2024/06/26 20:24

