

Event Sequence Delete

Version 6

Event Sequence Delete command will delete an event sequence. The command will only delete the series of events that were defined using the **Event Define Event Sequence** Command.

Example

Event_Define_Event_Sequence

```
/EVENT_SEQUENCE_NAME=RSTANCE  
/EVENT_SEQUENCE=RON+ROFF  
! /EXCLUDE_EVENTS=  
! /INSIDE_OF_SEQUENCE=  
/OFFSET_FROM_START=-1  
! /OFFSET_FROM_END=  
/OFFSET_BY=FRAMES  
/EVENT_SEQUENCE_INSTANCE=0
```

Event_Sequence_Delete

```
/EVENT_SEQUENCE_NAME=RSTANCE
```

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:event_commands:event_sequence_delete&rev=1721230935

Last update: 2024/07/17 15:42

