

Event Sequence Delete

Overview

Pipeline Command: Version 6

Event Sequence Delete command will delete an event sequence. The command will only delete the series of events that were defined using the **Event Define Event Sequence** Command.

```
Event_Sequence_Delete  
/EVENT_SEQUENCE_NAME=  
;
```

Example

```
Event_Define_Event_Sequence  
/EVENT_SEQUENCE_NAME=RSTANCE  
/EVENT_SEQUENCE=RON+ROFF  
! /EXCLUDE_EVENTS=  
! /INSIDE_OF_SEQUENCE=  
/OFFSET_FROM_START=-1  
! /OFFSET_FROM_END=  
/OFFSET_BY=FRAMES  
/EVENT_SEQUENCE_INSTANCE=0
```

```
Event_Sequence_Delete  
/EVENT_SEQUENCE_NAME=RSTANCE
```

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:event_commands:event_sequence_delete&rev=1741026064

Last update: 2025/03/03 18:21

