

Event Sequence Delete

Event Sequence Delete

Overview

The ``Event_Sequence_Delete`` command **removes a previously defined event sequence** from the motion file.

- Event sequences are often used to group event pairs (e.g., LHS to LTO) into meaningful motion phases. This command allows you to cleanly remove those sequences if no longer needed or if you're regenerating them.

Use cases:

- Cleaning up a trial before re-defining new event sequences.
- Removing incorrect or outdated gait cycle definitions.

Requirements:

- The event sequence to be deleted must already exist, and its name must match exactly.

Pipeline Command

```
Event_Sequence_Delete  
/EVENT_SEQUENCE_NAME=  
;
```

Command Parameters

Parameter	Description
<code>/EVENT_SEQUENCE_NAME=</code>	The name of the event sequence you want to delete (e.g., <code>`LeftGaitCycle`</code>).

Dialog

Due to the simplicity of this command, there is no dialog box and it can be edited within the text editor.

Examples

Simple Example

Delete the gait cycle sequence called `LeftStride`:

```
Event_Sequence_Delete
/EVENT_SEQUENCE_NAME=LeftStride
;
```

Complex Example

Use in a workflow where you are replacing event sequences:

```
Event_Sequence_Delete
/EVENT_SEQUENCE_NAME=RightCycle
;

Event_Define_Event_Sequence
/EVENT_SEQUENCE_NAME=RightCycle
/EVENT_SEQUENCE=RHS, RT0
/EXCLUDE_EVENTS=
/INSIDE_OF_SEQUENCE=
/OFFSET_FROM_START=
/OFFSET_FROM_END=
/OFFSET_BY=PERCENT
/EVENT_SEQUENCE_INSTANCE=
;
```

This removes the old `RightCycle` definition, then redefines it with a new event structure.

Notes

- Deleting an event sequence **does not delete the events themselves**, only the sequence structure. - Useful during iterative workflows where event phases are updated repeatedly. - Be cautious: if other commands depend on the sequence you remove, those commands may fail or behave unexpectedly.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:event_commands:event_sequence_delete&rev=1745857161

Last update: 2025/04/28 16:19

